

Distributed by



MASTERING IMPLANT REMOVAL

Be in control of any situation of
implant removal. OPERACE.
Flexible and efficient.



OPERACE
by PB Swiss Tools

 100% Made in Switzerland

This publication is not intended for distribution in the USA.

SURGICAL TECHNIQUE

PB SWISS TOOLS



TABLE OF CONTENTS

OPERACE	Innovation and Swiss quality	2
	Great choice for operating surgeons	3
APPLICATIONS	Applications	5
PRODUCT OVERVIEW	Sterile single-use inserts	7
	Instruments for multiple use	9
	OPERACE Instrument Set	10
SURGICAL TECHNIQUE	Using OPERACE correctly	13
	Procedure 1 Removal of intact screws that can be loosened by Hand	14
	Procedure 2 Removal of screws with a damaged screw recess	16
	Procedure 3 Removal of locking head screws stuck in the plate which cannot be removed by any of the above procedures	17
	Procedure 4a Removal of broken, projecting screws or screws processed according to procedure 3	18
	Procedure 4b Removal of buried broken screw shafts	19
	OPERACE Instrument Set MINI	20
	OPERACE Instrument Set SMALL	21
PRODUCT LINE	OPERACE Instrument Set LARGE	22
	Important Instructions General Instructions Instructions for single-use products Instructions for reprocessable products	23
IMPORTANT INSTRUCTIONS	Interpretation of Symbols	24

INNOVATION AND SWISS QUALITY

The OPERACE instrument set is used to remove implant screws made of titanium, titanium alloys and stainless steel – particularly damaged screws.

OPERACE is the result of the combined innovative power of PB Swiss Tools, the Swiss manufacturer of high-quality tools, and leading manufacturers of orthopedic implants and instruments.

The practical instrument set was developed in close cooperation with doctors and specialists. It offers a fast solution for the extraction of implant screws. The intelligent combination of durable handles and a range of single-use inserts allows for precise, safe and fast extractions of any type of screw.

OPERACE comes in three practical sizes:

- Instrument set MINI for mini fragment
- Instrument set SMALL for small fragment
- Instrument set LARGE for large fragment

A detailed product overview can be found on page 7 of this Surgical Technique.



GREAT CHOICE FOR SURGEONS

- **Convenience** With OPERACE you are equipped to remove standard implant screws from various manufacturers. Clearly structured, flexible and user-friendly through the intelligent color-coding of the three set sizes for implant screws with mini, small and large fragment.
- **Safety** Sterile single-use inserts.
- **Efficiency** The ergonomic PB Swiss Tools handles fit the hand conveniently and allow for direct power transmission.
- **Economic** Time-consuming and expensive resterilization of instruments that are not required is a thing of the past. Resulting in an efficient surgical procedure and clear cost transparency.



OPERACE
Instrument set SMALL
for small fragment









APPLICATIONS

For the removal of undamaged and damaged or broken screws, in particular

- Angular stable screws
- Cortex screws
- Cancellous bone screws
- Shaft screws
- Cannulated Screws
- Locking screws
- Locking bolts

made of titanium, titanium alloys and stainless steel with the following standard drives¹⁾:

Hex	Torx® / Stardrive®	Square/ Robertson	Cruciform	Slotted	Phillips
Screw head  HEX	Screw head  T	Screw head  SQR	Screw head  CR	Screw head  SL	Screw head  PH
Screw recess size 1.3 mm, 1.5 mm, 1.8 mm, 2.0 mm, 2.5 mm, 3.0 mm, 3.5 mm, 4.0 mm, 3/16", 5.0 mm	Screw recess size T4, T5, T6, T7, T8, T9, T10, T15, T20, T25, T30, T40	Screw recess size 0.7 mm, 1.0 mm, 1.2 mm, 1.5 mm, 1.8 mm, 2.2 mm, 2.3 mm	Screw recess size 0.3 mm, 0.4 mm, 0.5 mm, 0.6 mm, 0.8 mm, 1.0 mm	Screw recess size 0.3 mm, 0.4 mm, 0.5 mm, 0.6 mm, 0.8 mm, 1.0 mm	Screw recess size PH1

The system is organized into three color-coded sets based on screw diameter:

- Instrument set MINI for mini fragment, usually for screw diameters 0.9–2.0, yellow color-coding
- Instrument set SMALL for small fragment, usually for screw diameters 2.3–4.0, blue color-coding
- Instrument set LARGE for large fragment, usually for screw diameters 4.2–8.0, green color-coding

¹⁾ The screwdriver inserts in OPERACE correspond to standard orthopedic screwdriver inserts and are at least compatible with screws that satisfy the following specifications: ASTM F 543, ISO 5835, ISO 10664 and ISO 9268

HIGH-QUALITY SPECIAL STEEL

The drills are made from HSS (High Speed Steel). The special alloy of PB Swiss Tools is used for all other inserts. This alloy is particularly hard and tougher than other steels, with a higher breaking strength.



STERILE SINGLE-USE INSERTS

Every OPERACE instrument set contains single use inserts which can be used with OPERACE handles, cross-handles and extensions or power tools with adapters.

The single use inserts, extraction screws, drills and reamers are supplied in sterile form in double blister packaging and are intended for single use only.



OPERACE screwdriver inserts

Removal of intact screws with the following drives:

- Hexagonal socket (Hex)
- Torx®/Stardrive®
- Square socket (Robertson)
- Cruciform
- Slotted
- Phillips



OPERACE extraction screws

With left-hand thread, and available in a range of sizes, for the removal of screws with a damaged hexagonal socket, Torx®/Stardrive® or square socket drive.



OPERACE extraction drill bit

Removal of stuck locking-head screws. The screw head is separated from the screw shaft by drilling, enabling the plate to be removed with the separated screw head.



OPERACE extraction reamer


Left-turning, removal of drilled or broken screws. Three operations in one: reaming, gripping, extracting. Various sizes matched to their thread diameter.

Notes

- Sterile-packed, for single use. Inserts may not be reprocessed.
- After determining the correct insert size, remove the sterile insert from the packaging using aseptic surgical techniques.
- Do not resterilize and do not reuse if the packaging is damaged or torn.
- Unused components may not be resterilized.

ERGONOMIC HANDLES

Proven expertise gained in the product design of the industrial tools of PB Swiss Tools aided the development of the OPERACE handles and cross-handles. The result was two ergonomic handles for medical engineering.



INSTRUMENTS FOR MULTIPLE USE

The OPERACE handles, cross-handles, extensions and AO adapters are supplied in **non-sterile** form. They are intended for multiple use (reprocessing).

OPERACE handle

- For fast rotation and universal use
- The ergonomic groove shape rests comfortably in the hand and allows high torques to be transmitted
- The large rounded end of the handle enables considerable axial pressure to be applied

OPERACE cross-handle

- For the safe application of extremely high torques
- Optimally adjusted to fit the hand perfectly
- Thanks to the slant angle, the wrist remains straight and is protected when using the screwdriver with high force



OPERACE handles with quick coupling

The handles and cross-handles with quick coupling are designed to accept all OPERACE screwdriver inserts, extraction screws and extensions.



OPERACE extensions with quick coupling

If required, the extensions can be inserted between the handle and the insert to facilitate the removal of deeper implanted screws.



OPERACE couplings with AO/Synthes adapter for power tools

Couplings with an AO/Synthes adapter can be used as a link between power tools and all OPERACE screwdriver inserts, extraction screws and extensions.

Important

- The handles, extensions and couplings are intended for the loosening and unscrewing of implant screws.
- Blows and bending loads should be avoided since they can lead to instrument damage or breakage.
- The reprocessing must be carried out as described in the section on reprocessing.

OPERACE INSTRUMENT SET

Three set sizes with practical color-coding

The color-coding at the end of the handles corresponds to the color-coded size classification of the inserts (MINI, SMALL or LARGE).

The purpose of the size specific handles for the screwdriver inserts i.e. the sizes MINI, SMALL and LARGE, is to ensure that the achievable manual torque of the handle is matched to the torque appropriate to the screw size.


Instrument set MINI for mini fragment

Yellow color-coding
Screw diameters Ø 0.9–2.0 mm
Quick coupling with one ring



Instrument set SMALL for small fragment

Blue color-coding
Screw diameters Ø 2.3–4.0 mm
Quick coupling with two rings



Instrument set LARGE for large fragment


Green color-coding
Screw diameters Ø 4.2–8.0 mm
Quick coupling with three rings



Insert insertion/removal



Pull back quick coupling



Insert or remove insert



Release quick coupling
to lock in place

Note

- Check for correct engagement by pulling on the screwdriver insert

Storage container

The handy storage container provides space for the sterile inserts and processed (sterile) handles of one set. The containers can be optionally labeled MINI (yellow), SMALL (blue) or LARGE (green).



81006

SURGICAL TECHNIQUE



USING OPERACE CORRECTLY

Extracting screws requires experience and expertise. The OPERACE instrument set may be used only by trained operators.

The range of extraction products can be used for the following four procedures:

1. Removal of intact screws.
2. Removal of screws with a damaged screw recess.
3. Removal of locking head screws stuck in the plate.
4. Removal of broken screws and of screws processed according to procedure 3.



Animation

Important

- The correct instrument sizes to remove screws must be carefully chosen to facilitate screw extraction and to reduce the risk of further screw damage during the procedure.
- Supplementing the information in this manual, instruction in the surgical use of these instruments by a surgeon experienced in their handling is strongly recommended.

PROCEDURE 1

Removal of intact screws that
can be loosened by hand

Instruments used

- Handle or cross-handle of the appropriate size.
- Screwdriver insert, sterile-packed, single-use.

Option

- Extensions for the removal of deeply inserted screws
- Adapter with AO coupling for the use of instruments with power tools.

Note

- Start with the handle.
- If the screw cannot be rotated, use the cross-handle.

Procedure

1. Select screwdriver insert of the appropriate size and shape according to table 1.
2. Carefully clean screw recess.
3. Lock screwdriver insert in the handle coupling.
4. Insert the screwdriver insert into the screw recess.
5. **Loosen screw manually** by turning counterclockwise.

Option

- After it is loosened, remove the screw with the adapter with AO coupling and power tool

Important

- Use slot screwdriver inserts manually only
- Ensure that the screwdriver insert is fully inserted into the screw recess, otherwise the insert may spin in the recess.
















Table 1

Compatibility matrix based on screw diameter

					
0.9	0.9	0.9	0.9	0.9	0.9
1.0	1.0	1.0	1.0	1.0	1.0
1.2	1.2	1.2	1.2	1.2	1.2
1.3	1.3	1.3	1.3	1.3	1.3
1.4	1.4	1.4	1.4	1.4	1.4
1.5	1.5	1.5	1.5	1.5	1.5
1.6	1.6	1.6	1.6	1.6	1.6
1.7	1.7	1.7	1.7	1.7	1.7
1.8	1.8	1.8	1.8	1.8	1.8
1.9	1.9	1.9	1.9	1.9	1.9
2.0	2.0	2.0	2.0	2.0	2.0
2.1	2.1	2.1	2.1	2.1	2.1
2.3	2.3	2.3	2.3	2.3	2.3
2.4	2.4	2.4	2.4	2.4	2.4
2.5	2.5	2.5	2.5	2.5	2.5
2.6	2.6	2.6	2.6	2.6	2.6
2.7	2.7	2.7	2.7	2.7	2.7
2.8	2.8	2.8	2.8	2.8	2.8
3.0	3.0	3.0	3.0	3.0	3.0
3.2	3.2	3.2	3.2	3.2	3.2
3.3	3.3	3.3	3.3	3.3	3.3
3.4	3.4	3.4	3.4	3.4	3.4
3.5	3.5	3.5	3.5	3.5	3.5
3.7	3.7	3.7	3.7	3.7	3.7
3.8	3.8	3.8	3.8	3.8	3.8
3.9	3.9	3.9	3.9	3.9	3.9
4.0	4.0	4.0	4.0	4.0	4.0
4.2	4.2	4.2	4.2	4.2	4.2
4.3	4.3	4.3	4.3	4.3	4.3
4.5	4.5	4.5	4.5	4.5	4.5
4.6	4.6	4.6	4.6	4.6	4.6
4.7	4.7	4.7	4.7	4.7	4.7
4.8	4.8	4.8	4.8	4.8	4.8
4.9	4.9	4.9	4.9	4.9	4.9
5.0	5.0	5.0	5.0	5.0	5.0
5.3	5.3	5.3	5.3	5.3	5.3
5.5	5.5	5.5	5.5	5.5	5.5
5.7	5.7	5.7	5.7	5.7	5.7
6.0	6.0	6.0	6.0	6.0	6.0
6.4	6.4	6.4	6.4	6.4	6.4
6.5	6.5	6.5	6.5	6.5	6.5
6.8	6.8	6.8	6.8	6.8	6.8
7.0	7.0	7.0	7.0	7.0	7.0
7.3	7.3	7.3	7.3	7.3	7.3
7.5	7.5	7.5	7.5	7.5	7.5
8.0	8.0	8.0	8.0	8.0	8.0
Operace Ref	Operace Ref	Operace Ref	Operace Ref	Operace Ref	Operace Ref
HEX 1.3 80196	T4 80069	SQR 0.7 80036	CR 0.3 80055	SL 0.3 80061	PH 1 80086
HEX 1.5 80046	T5 80070	SQR 1.0 80037	CR 0.4 80056	SL 0.4 80062	
HEX 1.8 80051	T6 80071	SQR 1.2 80044	CR 0.5 80060	SL 0.5 80065	
HEX 2.0 80054	T7 80078	SQR 1.5 80194	CR 0.6 80057	SL 0.6 80063	
HEX 2.5 80048	T8 80073	SQR 1.8 80039	CR 0.8 80058	SL 0.8 80064	
HEX 3.0 80197	T9 80078	SQR 2.2 80195	CR 1.0 80059	SL 1.0 80066	
HEX 3.5 80050	T10 80078				
HEX 4.0 80052	T15 80075				
HEX 3/16" 80053	T20 80077				
HEX 5.0 80168	T25 80079				
	T30 80081				
	T40 80083				

- Most commonly used screw recess
- Other possible screw recesses

The table provides an overview of the screw recesses used for the respective screw diameters. Cells highlighted in black denote the screw recesses most commonly used in medicine with the correspondingly recommended screwdriver inserts. Other screw recesses are used depending on the manufacturer in each case. These and the correspondingly recommended screwdriver inserts are highlighted in gray.

	HEX 1.3, HEX 1.5, HEX 1.8, HEX 2.0, HEX 2.5, HEX 3.0, HEX 3.5, HEX 4.0, HEX 3/16" (4.7), HEX 5.0
	T4, T5, T6, T7, T8, T9, T10, T15, T20, T25, T30, T40
	SQR 0.7, SQR 1.0, SQR 1.2, SQR 1.5, SQR 1.8, SQR 2.2, SQR 2.3
	CR 0.3, CR 0.4, CR 0.5, CR 0.6, CR 0.8, CR 1.0
	SL 0.3, SL 0.4, SL 0.5, SL 0.6, SL 0.8, SL 1.0
	PH1

Note

- Recommendation for the selection of screwdriver inserts, with no guarantee of accuracy or completeness.

PROCEDURE 2

Removal of screws with a damaged screw recess

Instruments used

- Handle or cross-handle of the appropriate size.
- Extraction screw of the appropriate size, sterile-packed, single-use.

Procedure

1. Select extraction screw of the appropriate size matched to the screw recess size according to table 2.
2. Lock extraction screw in the handle coupling.
3. Start turning the extraction screw **to the left, counterclockwise**, in the same axis as the screw to be removed.
4. Continue turning, applying constant pressure, until the conical left-handed thread is securely seated in the damaged recess and until sufficient torque is applied to unscrew the screw.
5. Then unscrew the screw, turning to the left.

Important

- Turn to the left, counterclockwise.
- Use extraction screws only for removing screws with a damaged drive recess.
- Use extraction screws manually only.

Note

- If the extraction screw spins in the screw recess, an attempt can be made to drill the screw recess with the corresponding extraction drill bit to anchor the extraction screw more deeply.

Note

- For HEX 5.0 mm, T30 and T40 screw drives, a slightly deeper hole can be drilled in the recess with extraction drill bit 4.0.




Table 2

Screw recess size				OPERACE REF	
	mm				
1.5–3.2	Hex 1.5	T6, T7	SQR 1.0, 1.2, 1.5	80018	Extraction screw 1.6
2.0–4.3	Hex 1.8, 2.0	T8, T9	SQR 1.8	80170	Extraction screw 2.0
2.3–6.0	Hex 2.5, 3.0	T10, T15	SQR 1.8, 2.2	80020	Extraction screw 2.6
3.5–8.0	Hex 3.5, 4.0	T20, T25	SQR 2.3	80022	Extraction screw 3.5

PROCEDURE 3

Removal of locking head screws stuck in the plate which cannot be removed by any of the above procedures

Instruments used

- Power tool with Jacobs chuck.
- Extraction drill bit of the appropriate size, sterile-packed, single-use.

Procedure

1. Select extraction drill bit of the appropriate size according to table 3.
2. Position the drill bit in the screw recess and start turning the drill bit matching the screw diameter **to the right, clockwise**, in the same axis as the screw to be removed and applying **very little pressure**.
3. Only drill down until the underside of the implant plate is reached. This either loosens the screw head from the shaft or else weakens it sufficiently to cause it to break off when the plate is lifted.
4. Remove the plate.
5. Remove the screws according to procedure 4.

Important

- Turn to the right, clockwise.
- Apply the drill bit while it is already rotating, not while it is stationary.
- Apply very little pressure initially, only using the weight of the power tool, and start with a slow speed.
- Cool the bit and aspirate drill chips during the procedure. In order to avoid overheating, relieve the pressure on the drill from time to time.
- When the bit has reached the base of the recess, increase the pressure on the screw in order to obtain good chip formation.
- Do not use the bit to drill in bone or drill out broken screwdriver tips.





Table 3

Screw diameter in mm

Ø	OPERACE REF	
1.0–1.5	80023	Extraction drill bit 1.5
1.6–2.0	80123	Extraction drill bit 2.0
2.1–2.5	80024	Extraction drill bit 2.5
2.6–3.2	80025	Extraction drill bit 3.2
3.3–4.0	80019	Extraction drill bit 4.0
4.1–5.0	80026	Extraction drill bit 5.0
5.1–6.5	80027	Extraction drill bit 6.5
6.6–7.5	80021	Extraction drill bit 7.5

PROCEDURE 4A

Removal of broken, projecting screws or screws processed according to procedure 3

Instruments used

- Standard T-handle with Jacobs chuck or
- Power tool with Jacobs chuck.
- Extraction reamer of the appropriate size, sterile-packed, single-use.

Procedure

1. Select extraction reamer of the appropriate size matched to the screw recess according to table 4.
2. Locate the stationary extraction reamer and turn in the direction of the screw axis to ream over the screw shaft remaining in the bone. Apply **slight pressure** initially, **turning to the left, counterclockwise**.
3. As soon as the reamer grips the screw shaft, continue reaming with increased pressure until the conical left-handed thread is securely seated on the screw shaft.
4. **When unscrewing the reamer do not relieve the pressure, but maintain the constant axial pressure and direction of rotation.**

Note

- A manual procedure with a T-handle and Jacobs chuck is preferred.

Important

- Turn to the left, counterclockwise.
When unscrewing the reamer do not relieve the pressure, but maintain the constant axial pressure and direction of rotation.
- If a power tool is used, it should be kept at a very low speed.
- Cool the reamer and aspirate the drill chips during the procedure.






Table 4

Screw diameter in mm



Ø	OPERACE REF	
1.4–2.0	80009	Extraction reamer 2.0
2.1–2.5	80028	Extraction reamer 2.5
2.6–3.2	80029	Extraction reamer 3.2
3.3–4.0	80030	Extraction reamer 4.0
4.1–4.5	80033	Extraction reamer 4.5
4.6–5.8	80035	Extraction reamer 5.8
5.9–7.5	80015	Extraction reamer 7.5

PROCEDURE 4B

Removal of buried broken screw shafts

Instruments used

- Standard T-handle with Jacobs chuck or
- Power tool with Jacobs chuck.
- Extraction reamer of the appropriate size, sterile-packed, single-use.

Procedure



1. Select extraction reamer of the appropriate size matched to the screw recess according to table 4.
2. Locate the stationary extraction reamer perpendicular to the bone so that the teeth rest on the bone as evenly as possible. Apply **slight pressure** initially, **turning to the left, counterclockwise**.
3. As soon as the reamer grips, ream in the direction of the screw axis under image intensifier control until it has been guided through the screw shaft.
4. **Remove the reamer from the bone at regular intervals and free it of accumulated bone material.**
5. Continue reaming with increased pressure until the conical left-handed thread is securely seated on the screw shaft.
6. **When unscrewing the reamer do not relieve the pressure, but maintain the constant axial pressure and direction of rotation.**

Important

- Turn to the left, counterclockwise.
- When unscrewing the reamer do not relieve the pressure, but maintain the constant axial pressure and direction of rotation.
- If a power tool is used, it should be kept at a very low speed. If the screw grips, a manual procedure with a T-handle and Jacobs chuck is preferred.

















Note

Since the external diameter of *short-threaded cancellous bone screws* is smaller in the shaft than in the thread, smaller extraction reamers can be used if a break occurs in the shaft area. Although a reamer diameter smaller than that shown in the table 4 is used for this purpose, at least 10 mm of shaft without thread must be present.



OPERACE INSTRUMENT SET MINI



Reprocessable instruments		
	80001	Handle MINI
	80006	Extension MINI, for screwdriver inserts
Sterile single-use inserts		
		80196 Screwdriver insert for hexagonal socket screws, Hex 1.3
		80046 Screwdriver insert for hexagonal socket screws, Hex 1.5
		80069 Screwdriver insert for screws Torx® and Stardrive®, T4
		80070 Screwdriver insert for screws Torx® and Stardrive®, T5
		80071 Screwdriver insert for screws Torx® and Stardrive®, T6
		80055 Screwdriver insert for cruciform recess screws, Cr 0.3
		80056 Screwdriver insert for cruciform recess screws, Cr 0.4
		80060 Screwdriver insert for cruciform recess screws, Cr 0.5
		80036 Screwdriver insert for square socket screws, Sqr 0.7
		80037 Screwdriver insert for square socket screws, Sqr 1.0
		80061 Screwdriver insert for slotted screws, SI 0.3
		80062 Screwdriver insert for slotted screws, SI 0.4
		80065 Screwdriver insert for slotted screws, SI 0.5
		80018 Extraction screw, 1.6
		80023 Extraction drill bit, 1.5 f/screws Ø 1.0–1.5
		80123 Extraction drill bit, 2.0 f/screws Ø 1.6–2.0
		80009 Extraction reamer, 2.0 f/screws Ø 1.4–2.0
Storage		
	81006	Storage container, empty

OPERACE INSTRUMENT SET SMALL



Reprocessable instruments		
	80002	Handle SMALL
	80004	Cross handle SMALL
	80007	Extension SMALL, for screwdriver inserts
	80010	Pick-up coupling, SMALL, with AO/Synthes adapter
Sterile single-use inserts		
	80051	Screwdriver insert for hexagonal socket screws, Hex 1.8
	80054	Screwdriver insert for hexagonal socket screws, Hex 2.0
	80048	Screwdriver insert for hexagonal socket screws, Hex 2.5
	80178	Screwdriver insert for screws Torx® and Stardrive®, T7
	80073	Screwdriver insert for screws Torx® and Stardrive®, T8
	80198	Screwdriver insert for screws Torx® and Stardrive®, T9
	80078	Screwdriver insert for screws Torx® and Stardrive®, T10
	80075	Screwdriver insert for screws Torx® and Stardrive®, T15
	80077	Screwdriver insert for screws Torx® and Stardrive®, T20
	80057	Screwdriver insert for cruciform recess screws, Cr 0.6
	80058	Screwdriver insert for cruciform recess screws, Cr 0.8
	80059	Screwdriver insert for cruciform recess screws, Cr 1.0
	80086	Screwdriver insert for Phillips screw, Ph 1
	80044	Screwdriver insert for square socket screws, Sqr 1.2
	80194	Screwdriver insert for square socket screws, Sqr 1.5
	80039	Screwdriver insert for square socket screws, Sqr 1.8
	80195	Screwdriver insert for square socket screws, Sqr 2.2
	80041	Screwdriver insert for square socket screws, Sqr 2.3
	80063	Screwdriver insert for slotted screws, Sl 0.6
	80064	Screwdriver insert for slotted screws, Sl 0.8
	80066	Screwdriver insert for slotted screws, Sl 1.0
	80170	Extraction screw, 2.0
	80020	Extraction screw, 2.6
	80024	Extraction drill bit, 2.5 f/ screws Ø 2.1–2.5
	80025	Extraction drill bit, 3.2 f/ screws Ø 2.6–3.2
	80019	Extraction drill bit, 4.0 f/ screws Ø 3.3–4.0
	80028	Extraction reamer, 2.5 f/ screws Ø 2.1–2.5
	80029	Extraction reamer, 3.2 f/ screws Ø 2.6–3.2
	80030	Extraction reamer, 4.0 f/ screws Ø 3.3–4.0
Storage		
	81006	Storage container, empty

OPERACE INSTRUMENT SET LARGE



Reprocessable instruments

80003	Handle LARGE
80005	Cross handle LARGE
80008	Extension LARGE, for screwdriver inserts
80011	Pick-up coupling, LARGE, with AO/Synthes adapter

Sterile single-use inserts

		80197	Screwdriver insert for hexagonal socket screws, Hex 3.0
		80050	Screwdriver insert for hexagonal socket screws, Hex 3.5
		80052	Screwdriver insert for hexagonal socket screws, Hex 4.0
		80053	Screwdriver insert for hexagonal socket screws, Hex 3/16" (4.7)
		80168	Screwdriver insert for hexagonal socket screws, Hex 5.0
		80079	Screwdriver insert for screws Torx® and Stardrive®, T25
		80081	Screwdriver insert for screws Torx® and Stardrive®, T30
		80083	Screwdriver insert for screws Torx® and Stardrive®, T40
		80022	Extraction screw, 3.5
		80026	Extraction drill bit, 5.0 f/ screws Ø 4.1–5.0
		80027	Extraction drill bit, 6.5 f/ screws Ø 5.1–6.5
		80021	Extraction drill bit, 7.5 f/ screws Ø 6.6–7.5
		80033	Extraction reamer, 4.5 f/ screws Ø 4.1–4.5
		80035	Extraction reamer, 5.8 f/ screws Ø 4.6–5.8
		80015	Extraction reamer, 7.5 f/ screws Ø 5.9–7.5

Storage

81006	Storage container, empty
-------	--------------------------

IMPORTANT INSTRUCTIONS

GENERAL INSTRUCTIONS

Disposal

Used single-use inserts must be disposed of correctly.

Return

Surgical instruments or sterile inserts may be returned only after consulting the distributor. All returned products must be decontaminated or sealed in the unopened original packaging. Unopened products with an exceeded expiry date may not be returned.

Disclaimer

The recommendations for storage, care, maintenance, reprocessing and sterilization have been carefully checked, conform to international standard ISO 17664 and are considered to be appropriate. The person who actually reprocesses the instruments is responsible for achieving the desired results with the provided equipment, materials and personnel in the reprocessing facility. Please note the additional comprehensive information and instructions in the document: CLEANING AND STERILIZATION INSTRUCTIONS FOR REUSABLE INSTRUMENTS

The user is responsible for complications or other negative consequences which may result from reasons such as an incorrect indication or surgical technique, inappropriate material selection, inappropriate application or handling of the instruments, or any kind of application that is not described in the intended uses and any incorrect use, and cannot be blamed on the manufacturer, importer or supplier of PB Swiss Tools products. No liability is accepted in the event of failure to observe the instructions in this Surgical Technique.

Single-use products

Products intended for single use should never be reused. The reuse or reprocessing of single-use products can harm the product or jointly reprocessed products and/or lead to product failure. Furthermore, the reuse or reprocessing of single-use products increases the risk of contamination, potentially resulting in injury, illness or death of the patient or user.

INSTRUCTIONS FOR SINGLE-USE PRODUCTS

These instructions for use apply to all sterile inserts supplied by PB Swiss Tools GmbH and listed in this Surgical Technique.

- Screwdriver inserts
- Extraction Screws
- Extraction reamers
- Extraction drill bits

Sterile-packed, for single use

Inserts may not be reprocessed or resterilized.

After determining the correct insert size, remove the sterile insert from the packaging using aseptic surgical techniques. Do not resterilize and do not reuse if the packaging is damaged or torn.

Plasma coating

The color-coding of the inserts consists of a plasma coating. Slight deviations in color or shading result from the manufacturing process and do not affect the quality of the inserts.

INSTRUCTIONS FOR REPROCESSABLE PRODUCTS




























These instructions for use apply to all reprocessable surgical instruments supplied by PB Swiss Tools GmbH and listed in this Surgical Technique.

- Handles
- Cross-handles
- Extensions
- AO/Synthes couplings

The instruments are supplied in non-sterile form and must be reprocessed before every use. This also applies to the first use of reprocessable surgical instruments after delivery. Sterilization of the delivery packaging is not possible or permitted. Please note that only validated methods may be used for cleaning/disinfection and sterilization. The respective national and any applicable internal regulations must be observed. For further information on cleaning and sterilization see CLEANING AND STERILIZATION INSTRUCTIONS FOR REUSABLE INSTRUMENTS

www.pbswisstools.com/ifu

INTERPRETATION OF SYMBOLS

	Catalogue number
	Lot No./Batch code
	Serial number
	Manufacturer
	Date of manufacture
	Expiry date
	Non-sterile
	Sterilized by radiation
	For single use
	Do not resterilize
	Do not use if the sterile packaging is damaged
	Caution, consult the instruction for use
	Hex screw recess
	Torx®/Stardrive® screw recess
	Square/Robertson screw recess
	Cruciform screw recess
	Slotted screw recess
	Phillips screw recess
	Screw diameter
	Extraction screw
	Extraction drill bit
	Extraction reamer
	CE mark with identification number of the recognized inspection body
	Applying pressure
	Very little pressure
	Turn clockwise
	Turn counterclockwise
	Use manually only
	Cooling
	Slow speed

PB SWISS TOOLS

Legal Manufacturer:
PB Swiss Tools AG
Bahnhofstrasse 24
CH-3457 Wasen
info@pbswisstools.com

EC	REP
----	-----

MDSS GmbH
Schiffgraben 41
30175 Hannover
Germany



Information



Distributed by:
Synthes GmbH
Eimattstrasse 3
4436 Oberdorf
Switzerland
Tel: +41 61 965 61 11
Fax: +41 61 965 66 00
www.depuysynthes.com

This publication is not intended for distribution in the USA.

All surgical techniques are available as PDF files at
www.depuysynthes.com/ifu

CE 0123