


screwdriver insert compatibility matrix



- Most commonly used screw recess
- Other possible screw recesses

The table provides an overview of the screw recesses used for the respective screw diameters. Cells highlighted in black denote the screw recesses most commonly used in medicine with the correspondingly recommended screwdriver inserts. Other screw recesses are used depending on the manufacturer in each case. These and the correspondingly recommended screwdriver inserts are highlighted in gray.

- HEX 1.3, HEX 1.5, HEX 1.8, HEX 2.0, HEX 2.5, HEX 3.0, HEX 3.5, HEX 4.0, HEX 3/16" (4.7), HEX 5.0
- ⊕ T4, T5, T6, T7, T8, T9, T10, T15, T20, T25, T30, T40
- SQR 0.7, SQR 1.0, SQR 1.2, SQR 1.5, SQR 1.8, SQR 2.2, SQR 2.3
- ⊕ CR 0.3, CR 0.4, CR 0.5, CR 0.6, CR 0.8, CR 1.0
- SL 0.3, SL 0.4, SL 0.5, SL 0.6, SL 0.8, SL 1.0
- ⊕ PH1

Note

- Recommendation for the selection of screwdriver inserts, with no guarantee of accuracy or completeness.

Please note the comprehensive information and instructions in Surgical Technique.